

SNOOZA CHAIR

SUPPORTIVE RELAXATION CHAIR

Alternate seating that supports postural alignment and control

The innovative Snooza chair provides alternate seating that supports postural alignment and control when a child is out of a wheelchair. The Planar seat angles of the Snooza Chair use gravity to keep the user safe & posturally stable, and the high armrests help promote better breathing.

Standard configuration includes an adjustable headrest, central breathable upholstery, padded hip belt and a single wedge to alter the seat tilt angle.

Optional extras include padded tray, abductor pommel, padded mobile base, side infill pads (enabling seat width adjustment) and extra wedges to increase the tilt angle of the seat.



KEY FEATURES:

- Low centre of gravity offering stable seating position that is safe and comfortable
- Easy care, wipe-down medical grade vinyl with a breathable, washable centre wrap
- Carry handles that make transportation easy
- Large storage pockets located on both sides of the chair
- 4 playful colours to choose from (Apple Green, Bright Orange, Royal Blue & Hot Pink)


SNOOZA CHAIR - SUPPORTIVE RELAXATION CHAIR

OPTIONS:

- Pommel to control abduction
- Side infill pads to initially reduce seat width and allow for growth
- Extra wedges to increase the tilt angel of the seat
- Upholstered mobile base for convenience and manoeuvrability
- Leg rest to add comfort and aid stability of lower limbs



Snooza with Mobile Base



Gravity Aided Positioning



Snooza with Leg Rest



Apple Green



Bright Orange



Royal Blue



Hot Pink

TECHNICAL DATA:

Size	Seat Depth (cm)	Backrest Height (cm)	Seat Width at Rear (cm)	Seat Width at Front (cm)
1	26	60	20* - 25	27* - 32
2	32	70	23* - 28	30* - 34
3	40	75	27* - 32	37* - 42
4	45	80	30* - 34	40* - 44

*Width is adjustable when using the optional infill pads

