


Pride Zero Turn scooter specs.

Weight Capacity	136 kg (300 lbs.)
Maximum Speed ¹	Up to 6.4 km/h (4 mph)
Ground Clearance ²	3.3 cm (1.3") at motor
Turning Radius ²	97.2 cm (38.25")
Overall Length ²	116.2 cm (45.75")
Overall Width ²	50.8 (20")
Front Tyres	5 x 17.8 cm (2" x 7") solid
Rear Tyres	6.35 x 22.9 cm (2.5" x 9") solid
Range Per Charge ^{2,4}	Up to 16 km (10 miles)
Total Weight Without Batteries	54 kg (120 lbs.)
Heaviest Piece When Disassembled	21.8 kg (48 lbs.) (rear section)
Standard Seating	Type: Foldable, molded plastic Weight: 12 kg (27 lbs.) Material: Black vinyl Dimensions: Width: 43.2 cm (17") Depth: 43.2 cm (17")
Motor Output	24-volt DC motor, sealed mini-transaxle
Brakes	Electronic, regenerative and electromechanical
Battery Requirements ^{3,5}	(2) 12V deep-cycle Size: (2) 22 AH Weight: Battery box 15 kg (32 lbs.) 6 kg each
Battery Charger	Off-board, 3.5 amp
Delta Tiller	Yes
Colours	 Fire Opal Pearl White Sapphire Blue