


My chair is so fast
my friends have
to run to keep up.

— Wyatt



TiLite Pilot

Ready for take off



Specifications


Frame features	Carbon frame with 1" lower aluminium tubing	Seat back height	11" to 13" with standard back
Seat width	8" - 14"	Center of gravity adjustment	.5" - 3.5"
Seat depth	9" - 15"	Transit options	Transit Tie-Down, Transport
Front seat height	12" - 16.5"	Transport weight	4.9 kg*
Rear seat height	12" - 16"	Weight limit	43 kg transit, 54 kg standard
Front frame angle	Unlimited	Seat to footrest	3" - 14.5" without high mount parts

*Specifications may vary based upon configuration.
*Weight varies by configuration and does not include rear wheels.

Frame colours



Anodized colours



permobil

Permobil Australia | 1300 845 483 | info.au@permobil.com | permobil.com.au

Pilot Flyer AU Rev A

permobil



Optimal rear wheel access

Built in forward growth

Ready to grow

The TiLite Pilot has built-in forward growth adjustment to grow with your child and maintain optimal rear wheel access. This gives kids the most efficient propulsion possible without needing to adjust or relearn their push mechanics. The TiLite Pilot lets kids cruise their world like never before.

Designed to move

The clean frame design features a seat pan and side panels that provide easy access to the adjustment bolts, making changes simple and intuitive. Unlike on other chairs, all of these adjustments can be done without compromising access to the rear wheels.

Easy to use wheel locks are within reach



Frame designed for maximum stability

Lightweight for easier transport

Cleared for landing

The ultra-lightweight design and carbon fiber frame panels make it 50% lighter and a breeze for kids to navigate all the adventures of their great big world. Quick release wheels and transit tie downs mean that bringing the pilot everywhere is also simple for mum and dad.

