

TiLite Pilot Specifications

Specifications

Frame features	Carbon frame with 1" lower aluminium tubing	Seat back height	11" to 13" <i>with standard back</i>
Seat width	8" - 14"	Center of gravity adjustment	.5" - 3.5"
Seat depth	9" - 15"	Transit options	Transit Tie-Down, Transport
Front seat height	12" - 16.5"	Transport weight	4.9 kg*
Rear seat height	12" - 16"	Weight limit	43 kg transit, 54 kg standard
Front frame angle	Unlimited	Seat to footrest	3" - 14.5" <i>without high mount parts</i>

*Specifications may vary based upon configuration.
*Weight varies by configuration and does not include rear wheels.